




Bariatric Patient Handling System

Purchase and Rental packages

The Polaris Medical free-standing gantry Bariatric System with its unique anti-tilting device* provides the complete solution to moving and handling patients up to 980 lbs!

Available to institutions on a purchase or rental basis, offering a flexible cost effective solution to manual handling needs.



* Patent Pending






Polaris Medical Inc.

www.polarismedicalinc.com

Option 1

Dual Hoist System & X-Y Gantry


(Lifting capacity 445kgs (70 stones))



Option 2

Single Hoist System & X-Y Gantry

(Lifting capacity 272kgs (42 stones))



Purchase & Rental Packages

Standard package includes:

Installation & Commissioning of system

- Survey prior to installation
- 24 hour technical service and support
- Modular Gantry sized to suit location
- 2 Standard Slings
Small, medium, large, extra large options - to be specified at time of order.

Rental package also includes:

- Removal on completion of rental
- Installation of system within 48 hours from receipt of order.

Optional for rental or purchase

- Weight Monitor for patient monitoring
- Physiotherapy and rehabilitation requirements can be accommodated subject to survey.
- Bariatric Blanket Sling (Order code SLBB) Maximum Load 445kg

Distributed by:
Polaris Medical Group

T 905.467.5152
F 905.333.0407
TF 877-377-9988

info@polarismedicalinc.com
www.polarismedicalinc.com

