



Patent Pending



...An Uplifting Experience!

The New HoverJack™ Air Patient Lift from HoverTech International allows a patient to be lifted from the floor to bed or stretcher level in a supine position, eliminating discomfort and the risk of injury to patient and staff.

The patient can be transported on the inflated HoverJack™ to the bedside with an optional friction-reducing material attached to the underside of the lift.

Once the HoverJack™ is inflated, the same air supply unit can be used with the HoverMatt® transfer mattress to safely transfer the patient to a bed or stretcher.



HoverTech International
800-471-2776
hovermatt@earthlink.net
www.hovermatt.com