



Polaris Medical Inc.

www.polarismedicalinc.com

The Curved Transfer Board is a very hard wearing polyethylene transfer board for seated independent or assisted transfers.

The board is shaped to facilitate positioning and features anti-slip pads on the underside for safer location and a convenient carrying handle.

The board has a max. user weight of 120Kgs and can be easily cleaned with soapy water/alcohol/disinfectant and a non-abrasive cloth.



- High Quality
- Easy to clean
- Unbreakable
- Chip resistant
- Made from recyclable Polyethylene



Polaris Medical Inc.
Burlington, Ontario
1-877-377-9988 T
1-905-333-0407 F
info@polarismedicalinc.com

1-905-635-7428